

SANITAS[®]

Synergy Vet



BROAD SPECTRUM DISINFECTANT WITHOUT ALDEHYDES



BACTERIA



VIRUSES



FUNGI

Latest generation disinfectant that combines two active ingredients that together exert a broad spectrum disinfectant action with reduced in-use doses: its excellent compatibility with all kinds of materials makes it a versatile product for the complete disinfection of livestock facilities and premises, as well as to disinfect all kinds of livestock transport vehicles, veterinary clinics, (cages, etc) and slaughterhouses (animal stabling area). It is also authorised for use in the food industry.



SANITAS®

Synergy Vet


DISINFECTANTS



COMPOSITION

Glutaraldehyde	8.5%
Didecyl dimethyl ammonium chloride	11.9%
Glycolic acid	0.7%
Excipient qs.	100%

FORMATS AVAILABLE



YouTube

To request the technical specifications and safety sheet, see demonstration videos or for further information scan this code:



INDICATIONS, RECOMMENDED DOSES AND USE:

SANITAS® SYNERGY VET can be used to disinfect all kinds of surfaces, whether by brushing or spraying since the efficacy is identical. It can also be used in drive-through disinfection baths.

The product will be applied at the recommended dose and left to dry on the application site or if this is not possible, leaving it to act for at least 30 minutes. Utensils or machinery can be disinfected by immersion in a solution at the indicated dose of the product for a minimum of half an hour and subsequently rinsed with large amounts of water if its size so permits, or alternatively disinfected by spraying.

Recommended doses:

Dilute the product in cold water in the following proportions by using the accompanying measuring cap:

Premises and buildings

Microorganism	Condition	Time (min)	Dosage
Bacteria	Dirty	30	0.1%
Viruses	Dirty	5	0.75%
Fungi and yeasts	Clean	30	1%
Aerial disinfection by thermal fogging against bacteria, fungi and viruses	Clean	30	1%
	Dirty	30	25% (1 l of SANITAS SYNERGY VET per every 1000 m ³)
Disinfection against ASF	Dirty	5	0.2%
Disinfection against avian influenza	Clean	5	0.2%

Utensils: 1:200 (50 ml SANITAS® SYNERGY VET in 10 L of water).

Means of transport: 1:200 (50 ml SANITAS® SYNERGY VET in 10 L of water).

Drive-through disinfection baths: 1:100 (100ml of SANITAS® SYNERGY VET in 10L of water)

Apply in the following proportions:

100 ml dilution per square meter on non-porous surfaces.

300 ml dilution per square meter on porous surfaces.

PHYSICO-CHEMICAL CHARACTERISTICS:

Formulation: Concentrated liquid.

Appearance: Clear intense blue liquid with a characteristic odor.

Density (20° C): 1.010-1.025 g/cc.

Solubility in water: soluble.

pH: 4,5 (concentrated product).

SHELF LIFE:

SANITAS® SYNERGY VET has a shelf life of 4 years.

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

Laboratorios Zotal S.L.U.

Ctra. Nacional 630, km. 809.

41900 Camas (Sevilla) SPAIN

Tel.:+34 954 390 204 Fax:+34 954 395 516

Email: zotal@zotal.com



efr

USER INFORMATION: +34 954 395 297

www.zotal.com